



# PLASTIC SURGERY INSTRUMENTS



**M.F. ELAHI MEDICAL CO**  
SERVING WORLDWIDE SURGICAL INSTRUMENTS

# TECHNICAL SUPPORT & WARRANTY INFORMATION



M.F. ELAHI MEDICAL CO

## TECHNICAL SERVICE

At MF Elahi Medical Co., precision and performance are at the core of every instrument we deliver.

We provide full technical assistance to ensure:

- Precision manufacturing for microsurgical procedures
- Product guidance and instructions for safe handling
- Instrument care tips to extend lifespan
- Compatibility with international sterilization standards

## AFTER-SALES SUPPORT

We are committed to a seamless customer experience:

- Fast response to product-related inquiries
- Support with handling, storage, and usage protocols
- Product documentation & compliance certificates available upon request
- Dedicated contact via WhatsApp and Email

## WARRANTY POLICY

We offer a 1-year limited warranty on all Plastic Surgery Instruments:

- Covers manufacturing defects only
- Warranty valid only when instruments are used for their intended surgical purpose
- Misuse, modification, or improper sterilization voids warranty

## CERTIFICATIONS

MF Elahi Medical Co. products are manufactured under rigorous international standards and hold the following certifications:

**Certification Standard**

- ✓ CE Mark European Conformity for Medical Devices
- ✓ ISO 13485 Medical Devices Quality Management
- ✓ ISO 9001 General Quality Management
- ✓ GMP Compliant Good Manufacturing Practices

 **Need Assistance?**

Reach us any time:

 [zakriamfe@gmail.com](mailto:zakriamfe@gmail.com) / [mfe@inter.net.pk](mailto:mfe@inter.net.pk)

 [www.mfelahi.com](http://www.mfelahi.com)



PLASTIC SURGERY INSTRUMENTS	Page No.
FACE LIFT SCISSOR	2
DISSECTING SCISSOR	3
SCISSORS	4
TENOTOMY SCISSORS	5
EASY-CUT SCISSOR	6-7
DISSECTING FORCEPS	8-9
FACE LIFT FORCEPS	10
NEEDLE HOLDER	11-12
RETRACTORS & SKIN HOOKS	13
RETRACTORS	14
SELF RETAINING RETRACTORS	15
DISSECTORS	16
SCALPEL HANDLES & MARKING PENS	17
SKIN GRAFT KNIVES	18
AREOLA MARKERS	19
FIBER OPTIC BREAST RETRACTORS & LEADS	20



NASAL INSTRUMENTS	Page No.
NASAL SCISSORS	21-22
SEPTUM KNIFE	23
SEPTUM ELEVATOR	24
NASAL RETRACTORS	25-26
NASAL SPECULUM	27
NASAL SAWS, KNIVES & SEPTUM FORCEPS	28-29
MALLETS & GOUGES	30
NASAL OSTEOTOMES & CHISELS	31-33
NASAL RASP	34-36



ENDOSCOPIC BROWLIFT INSTRUMENTS	Page No.
BROWLIFT ELEVATORS	37-42



MAXILLO FACIAL INSTRUMENTS	Page No.
OSTEOTOMES	43-44
DISIMPACTION FORCEPS	45-46
SMITH RAMUS SEPARATOR	47
FRAGMENT FORCEPS	48



MICRO INSTRUMENTS	Page No.
VESSEL DILATOR	49-51
MICRO FORCEPS	52-53
MICRO SCISSORS	54
MICRO NEEDLE HOLDERS	55
ULTRA FINE MICRO INSTRUMENTS	56-59
INSTRUMENTS CONTAINER	60



- Tungsten carbide edges

Straight	Curved	Length
MFE-691-17	MFE-692-17	17cm/6 3/4"
MFE-691-19	MFE-692-19	19cm/7 1/2"
MFE-691-23	MFE-692-23	23cm/9"

Flat tips



**Gorney Face Lift Scissors**

- Tungsten carbide edges
- Open shanks

Straight	Curved	Length
MFE-693-17	MFE-694-17	17cm/6 3/4"
MFE-693-19	MFE-694-19	19cm/7 1/2"
MFE-693-23	MFE-694-23	23cm/9"

Flat tips



**Gorney-Freeman Face Lift Scissors**

- Tungsten carbide edges

- Curved

Code	Length
MFE-695-17	17cm/6 3/4"
MFE-695-19	19cm/7 1/2"
MFE-695-23	23cm/9"



**Kaye Face Lift Scissors**

- Tungsten carbide edges

- Curved

Code	Length
MFE-696-14	14cm/5 1/4"



**Aufrecht Scissors**

- Tungsten carbide edges

- Curved

- Tapered blades to ne tips

Code	Length
MFE-697-13	13cm/5"



**Goldman-Fox Scissors**



- Tungsten carbide edges
- Curved at tips

**Code**      **Length**

<b>MFE-698-13</b>	<b>13cm/5 1/4"</b>
<b>MFE-698-15</b>	<b>15cm/6"</b>
<b>MFE-698-18</b>	<b>18cm/7"</b>




**Ragnell (Kilner) Scissors**

- Tungsten carbide edges
- Curved Tepered at tips

**Code**      **Length**

<b>MFE-699-11</b>	<b>11.5cm/4 1/4"</b>
-------------------	----------------------



**Greenberg PAR Scissors**

- Hole in the blades for the suture to be threaded through
- Curved

**Code**      **Length**

<b>MFE-700-12</b>	<b>12cm/4 3/4"</b>
-------------------	--------------------



**Littler Plastic Surgery Scissors**

- Super Cut edges
- Curved

**Code**      **Length**

<b>MFE-701-14</b>	<b>14.5cm/5 3/4"</b>
-------------------	----------------------



**Kahn Face Lift Scissors**

- Super Cut edges
- Curved, round blunt tips

**Code**      **Length**

<b>MFE-702-13</b>	<b>13.5cm/5 1/4"</b>
-------------------	----------------------




**Fomon Saber Face Lift Scissors**



- Micro serrated cutting edge
- Straight

**Code**      **Length**

**MFE-703-11**      11.5cm/4½”



**Iris Scissors**

- Micro serrated cutting edge
- Curved

**Code**      **Length**

**MFE-704-11**      11.5cm/4½”



**Iris Scissors**

- Micro serrated cutting edge
- Straight

**Code**      **Length**

**MFE-705-11**      11.5cm/4½”



**Strabismus Scissors**

- Micro serrated cutting edge
- Curved

**Code**      **Length**

**MFE-706-11**      11.5cm/4½”



**Strabismus Scissors**

- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-707-12**      12.5cm/5”





**Waldram Scissors (with larger bows)**



- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-708-11**      **11.5cm/4½"**




**Stevens Tenotomy Scissors**

- Curved
- Micro serrated cutting edge
- Larger oversized rings

**Code**      **Length**

**MFE-709-12**      **12.5cm/5"**



**Stevens-Stille Pattern Tenotomy Scissors**

- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-710-13**      **13.5cm/5¼"**




**Jameson-Werber Scissors**

- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-711-15**      **15cm/6"**



**Jameson Scissors**

- Curved Strong pattern
- Micro serrated cutting edge

**Code**      **Length**

**MFE-712-15**      **15cm/6"**



**Jameson Scissors**



- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-713-14**      **14cm/5½"**



**Slim Baby Metzenbaum Easy-Cut Scissors**

- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-714-18**      **18cm/7"**



**Slim Metzenbaum Easy-Cut Scissors**

- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-715-14**      **4cm/5½"**



**Baby Metzenbaum Easy-Cut Scissors**

- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-716-18**      **18cm/7"**



**McIndoe / Metzenbaum Easy-Cut Scissors**

- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-717-18**      **18cm/7"**




**Crayford-Reynolds Scissors**



- Straight
- Micro serrated cutting edge

**Code**      **Length**

**MFE-718-14**      **14cm/5½"**



**Mayo Easy-Cut Scissors**

- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-719-14**      **14cm/5½"**




**Mayo Easy-Cut Scissors**

- Straight
- Micro serrated cutting edge

**Code**      **Length**

**MFE-720-17**      **17cm/6¾"**



**Mayo Easy-Cut Scissors**

- Curved
- Micro serrated cutting edge

**Code**      **Length**

**MFE-721-17**      **17cm/6¾"**



**Mayo Easy-Cut Scissors**




• **Toothed**  
**Code**

MFE-722-15

**Length**

15cm/6"



**Gillies Dissecting Forceps**

• **Non Toothed**  
**Code**

MFE-723-15

**Length**

15cm/6"



**McIndoes Dissecting Forceps**


• **1.5mm tips**

**Code**

MFE-724-15

**Length**

15cm/6"



**Debakey-Diethrich Tissue Forceps**


• **2mm tips**

**Code**

MFE-725-15

**Length**

15cm/6"



**Debakey Tissue Forceps**

• **1mm tips**

**Code**

MFE-726-12

**Length**

12.5cm/5"



**Bipolar Forceps**


• **1.5mm tips**

**Code**

MFE-727-15

**Length**

15cm/6"



**Bipolar Forceps**

• **Block end/Banana plug thing to machine end**

**Code**

MFE-728-03

**Length**

3mm





**Bipolar Cable**



• **Non Toothed**

Code	Length
MFE-729-12	12cm/4 <sup>3</sup> / <sub>4</sub> "



**Adson Forceps**

• **Toothed**

Code	Length
MFE-730-12	12cm/4 <sup>3</sup> / <sub>4</sub> "



**Adson Forceps**

• **Non Toothed**

Code	Length
MFE-731-12	12cm/4 <sup>3</sup> / <sub>4</sub> "



**Micro Adsons Forceps**

• **Toothed**

Code	Length
MFE-732-12	12cm/4 <sup>3</sup> / <sub>4</sub> "



**Micro Adsons Forceps**

• **7x7 Teeth**

Code	Length
MFE-733-12	12cm/4 <sup>3</sup> / <sub>4</sub> "



**Adson Brown Forceps**

• **2mm tips**

• **Atraumatic tips**

Code	Length
MFE-734-12	12cm/4 <sup>3</sup> / <sub>4</sub> "



**Adson-DeBakey Forceps**



- Diamond dusted jaws

Code	Length
MFE-735-12	12.5cm/5"



**Assumpcao Face Lift Forceps**

- Tungsten carbide jaws

Code	Length
MFE-736-12	12cm/4 <sup>3</sup> / <sub>4</sub> "



**Pitanguy Forceps**

- Tungsten carbide jaws

Code	Length
MFE-737-13	13cm/5"



**Marten Flap Forceps**

- Tungsten carbide jaws

Code	Length
MFE-738-15	15cm/6"
MFE-738-18	18cm/7"
MFE-738-20	20cm/8"



**Martel Flap Grasping Forceps**



- Diamond dusted jaws

Code	Length
MFE-739-12	12.5cm/5"



**Micron Fine Needle Holder**

- Tungsten carbide jaws

Code	Length
MFE-740-12	12cm/4 <sup>3</sup> / <sub>4</sub> "



**Derf Needle Holder**

- Tungsten carbide jaws

Code	Length
MFE-741-13	13cm/5"



**Halsey Needle Holder**

- Tungsten carbide jaws

Code	Length
MFE-742-15	15cm/6"
MFE-742-18	18cm/7"
MFE-742-20	20cm/8"




**Crile Wood Needle Holder**



- Tungsten carbide jaws


Code	Length
MFE-743-13	13cm/5¼"



**Nieverts Needle Holder**  
(o set ring)

- Tungsten carbide jaws

Code	Length
MFE-744-13	13cm/5"



**Foster Gillies Combined Needle Holder/Scissors**

- Tungsten carbide jaws

Code	Length
MFE-745-16	16cm/6¼"



**Gillies Combined Needle Holder/Scissors**



**Code Length**

**MFE-746-15 15cm/6"**



**Kilner Detachable Skin Hook**

**Code Length**

**MFE-746-18-S 18cm/7"-S**

**MFE-746-18-L 18cm/7"-L**



**Gillies Skin Hook**

**Code Length**

**MFE-747-16-S 16cm/6¼"-S**

**MFE-747-16-L 16cm/6¼"-L**



**Klienert-Kutz Skin Hook**

**Code Length**

**MFE-748-16 16cm/6¼"**



**Joseph Double Skin Hook**

- Semi-sharp prongs

**Code Length**

**MFE-749-13 13.5cm/5"**



**Rollet/Blair Mini-Catspaw Retractor**

- Semi-sharp prongs

**Code Length**

**MFE-750-16 16cm/6¼"**



**Catspaw Retractors**



**Code**      **Length**

MFE-751-15      15cm/6"



**Kilner-Miniature Pattern Catspaw Retractor**

**Code**      **Length**

MFE-752-18-S      18cm/7"-S

MFE-752-18-L      18cm/7"-L



**Senn-Miller Double Ended Retractor**

**Code**      **Length**

MFE-753-16-S      16cm/6¼"-S

MFE-753-16-L      16cm/6¼"-L



**Ragnell-Davies Retractor**

**Code**      **Length**

MFE-754-16      16cm/6¼"



**Mini-Langenbeck Retractor**

- Semi-sharp prongs

**Code**      **Length**

MFE-755-13      13.5cm/5"



**Desmarres Saddle Retractor**

- Semi-sharp prongs

**Code**      **Length**

MFE-756-16      16cm/6¼"




**Langenbeck Retractor**



- 4 x 4 teeth blunt
- 4mm depth prongs

**Code**      **Length**

**MFE-757-07**      **7cm/2<sup>3</sup>/<sub>4</sub>"**



**Alms Retractor**

- 2 x 3 blunt prongs

**Code**      **Length**

**MFE-758-11**      **11.5cm/4<sup>1</sup>/<sub>2</sub>"**



**Baby Weitlaner S/R Retractor**

- 3 x 4 sharp prongs

**Code**      **Length**

**MFE-759-14**      **14cm/5<sup>1</sup>/<sub>2</sub>"**



**Wests S/R Retractor**

- 2 x 2 sharp prongs
- 12mm depth Curved

**Code**      **Length**

**MFE-760-13**      **13cm/5<sup>1</sup>/<sub>4</sub>"**



**Mollison S/R Retractor - Child**

- 4 x 4 sharp prongs
- 14mm depth Curved

**Code**      **Length**

**MFE-761-17**      **17cm/6<sup>3</sup>/<sub>4</sub>"**



**Mollison S/R Retractor - Adult**



**Code**      **Length**

**MFE-762-18**      **18.5cm/7¼"**



**Freers Elevator/Dissector**

**Code**      **Length**

**MFE-763-18**      **18cm/7"**

**MFE-763-23**      **23cm/9"**

**MFE-763-25**      **25.5cm/10"**



**Watson Cheyne Dissector**

**Code**      **Length**

**MFE-764-19**      **19cm/7½"**

**MFE-764-23**      **23cm/9"**



**MacDonald Elevator/Dissector**

**Code**      **Length**

**MFE-765-21**      **21.5cm/8½"**



**Howarth Elevator**

**Code**      **Length**

**MFE-766-16**      **16cm/6¼"**



**Mitchels Trimmer**



**Code Length**

**MFE-767-13 13cm/5¼"**



**Barrons Scalpel Blade Handle No.3**

**Code Length**

**MFE-768-13 13cm/5¼"**



**Beaver Blade Holder**

**Code Length**

**MFE-769-14 14.5cm/5¾"**



**Ecko Mapping Pen**

**Code Length**

**MFE-770-15 15cm/6"**



**Sommerlad Marking Pen**

- Measuring range 0 to 25mm

**Code Length**

**MFE-771-07 7cm/2¾"**



**Castroviejo Measuring Caliper**

**Code**

**MFE-772**



**Ink Pot Stainless Steel**

**Code Length**

**MFE-773-15 15cm/6"**

**MFE-773-30 30cm/12"**



**Stainless Steel Measuring Ruler**




Right Handed	Left Handed	Length
MFE-774-30	MFE-775-30	30.5cm/12"



The roller bar remains smooth and stationary during use. Graft thickness is adjusted by loosening the locking nut, turning the adjuster nut to the desired setting, and then securing it by tightening the locking nut. The roller bar adjusts only in parallel, ensuring consistent graft thickness.

**Watson Skin Graft Knife**


Right Handed	Left Handed	Length
MFE-776-30	MFE-777-30	30.5cm/12"



The grooved roller bar both oscillates laterally and rotates. Graft thickness is adjusted by turning the adjuster nuts at the knife's end, which raises or lowers the roller bar.

**Braithwaite Skin Graft Knife**


Right Handed	Left Handed	Length
MFE-778-30	MFE-779-30	30.5cm/12"



The grooved roller bar oscillates laterally without rotating. Graft thickness adjustment follows the same method as the Braithwaite knife.

**Cobbett Skin Graft Knife**

Right Handed	Left Handed	Length
MFE-780-36	MFE-781-36	36.5cm/14 1/4"



The grooved roller bar both rotates and oscillates laterally. Graft thickness is adjusted using the thumb screws on top of the support pillars.

**Humby Skin Graft Knife**

Right Handed	Left Handed	Length
MFE-782-19	MFE-783-19	19cm/7 1/2"



**Silver Skin Graft Knife**

- 200mm x 100mm/8 x 4" Bevelled edge

Knife	Blade Code	Code	Size	Right Handed	Left Handed	Length
Cobbett	9940	MFE-784	200mm x 100mm/8 x 4"	MFE-785-19	MFE-786-19	19cm/7 1/2"
Watson	9941					
Braithwaite	9942					
Silver	9943					



**Tufnol Skin Graft Board**



Code	Diameter
MFE-787-36	36mm Ø
MFE-787-38	38mm Ø
MFE-787-40	40mm Ø
MFE-787-42	42mm Ø
MFE-787-45	45mm Ø



**Freeman Areola Marker**

Code	Diameter
MFE-788-38	38/42mm Ø
MFE-788-42	42/45mm Ø



**Freeman Double Ended Areola Marker**

Code	Diameter
MFE-789-40	40mm Ø
MFE-789-44	44mm Ø
MFE-789-48	48mm Ø
MFE-789-52	52mm Ø



**Grossman Areola Marker**


Code	Diameter
MFE-790-40	40mm Ø standard
MFE-790-50	50mm Ø large



**Galvao Nipple Marker**

- Accurately measures the amount of tissue designated for removal.
- Designed for use with the Freeman Marker.
- Featuring arms graduated in centimeters to facilitate precise marking of the tissue to be excised.


Code	Diameter
MFE-791-38	38mm Ø
MFE-791-42	42mm Ø
MFE-791-45	45mm Ø



**McKissock Keyhole Breast Reduction Marker**



- Tebbetts Fibre Optic Breast Retractor With fibre optic illumination and suction



- Fibre optic illumination with suction for smoke evacuation
- Angled Retractor tip with blunt teeth
- Polished mirror Finish on blade to reflect light inside the pocket

**With Teeth**

Code	Length	Width
MFE-792-24	90mm	24mm
MFE-792-30	90mm	30mm
MFE-792-36	150mm	36mm
MFE-792-190	190mm	36mm



Code	Blunt Length	Width
MFE-793-24	90mm	24mm
MFE-793-30	90mm	30mm
MFE-793-36	150mm	36mm
MFE-793-190	190mm	36mm



- Tungsten Carbide edges
- Double angled blades

**Code**      **Length**

**MFE-794-13**      **13.5cm/5¼"**



**Giunta Nasal Scissors**

- 40° angled blades

**Code**      **Length**

**MFE-795-10**      **10cm/4"**



**Daniels Nasal Scissors**

- Curved

**Code**      **Length**

**MFE-796-13**      **13.5cm/5¼"**



**Fomon Upper Lateral Scissors**

- Tungsten Carbide edges

**Code**      **Length**

**MFE-797-13**      **13cm/5"**



**Fomon Dorsal Scissors**

- Tungsten Carbide edges

**Code**      **Length**

**MFE-798-16**      **16cm/6¼"**



**Cottle Dorsal Scissors**



- Tungsten Carbide edges

Code	Length
MFE-799-18	18cm/7"



**Knight Nasal Scissors**

Code	Length
MFE-800-18	18cm/7"



**Heymann Nasal Scissors**

- Double action angled

Code	Length
MFE-801-20	20cm/8"



**Caplan Nasal Shears**



**Code Length**

**MFE-802-14 14cm/5½"**




**Cottle Septum Knife**

• **Straight**

**Code Length**

**MFE-803-14 14cm/5½"**



**Joseph button end Knife**

• **Curved**

**Code Length**


**MFE-804-16 16cm/6¾"**



**Convex button end Knife**

**Code Length**


**MFE-805-15 15cm/6"**



**Freer Mucosa Knife**

**Code Length**

**MFE-806-17 17.5cm/7"**



**Joseph Angled button end Knife**



- Double ended graduated

Code	Length
MFE-807-20	20cm/8"



**Code**      **Length**

MFE-808-20	20cm/8"
------------	---------



- Double ended with curette

**Code**      **Length**

MFE-809-21	21cm/8¼"
------------	----------



- Double ended

**Code**      **Length**

MFE-810-19	19cm/7½"
------------	----------



- Self retaining ratchet
- Adjustable to hold breast during surgery

**Code**

MFE-811
---------




**Mammoplast Breast Retractor**



- 85mm length

Code	Blunt Prongs	Code	Sharp Prongs
MFE-812-10	10mm	MFE-813-10	10mm
MFE-812-13	13mm	MFE-813-13	13mm
MFE-812-15	15mm	MFE-813-15	15mm



**Kilner Alae Retractor**

- 2 prong blunt ball tips
- Knife guard

Code	Length
MFE-814-20	20cm/8"



**Cottle-Neivert Alae Retractor**

- 2 prong blunt ball tips

Code	Length
MFE-815-16	16cm/6¼"



**Fomon Alae Retractor**


- 2 prong blunt ball and sharp tips

Code	Length	
MFE-816-01	14cm/5½"	Fig.1
MFE-816-02	14cm/5½"	Fig.2
MFE-816-03	14cm/5½"	Fig.3




**Cottle Hook**

Fig.1



**Sharp Prongs**

Fig.2



**Ball on Right**

Fig.3



**Ball on Left**



**Code**

**Length**

**MFE-817-15**

**15cm/6"**



**Cottle Retractor**

- 45mm blade

**Code**

**Length**

**MFE-818-16**

**16cm/6<sup>1</sup>/<sub>4</sub>"**



**Aufrecht Nasal Retractor**

- Fibre Optic illumination
- 45mm blade

**Code**

**Length**

**MFE-819-16**

**16cm/6<sup>1</sup>/<sub>4</sub>"**



**Aufrecht Nasal Retractor**



- 13cm/5” Length
- Plain Blade

**Code**

**Blade**

**MFE-820-35**

**35mm**

**MFE-820-55**

**55mm**

**MFE-820-65**


**65mm**

**MFE-820-75**

**75mm**

**MFE-820-90**

**90mm**



**Killian Speculum**

- 14cm/5½” length
- Slender blade

**Code**

**Blade**

**MFE-821-20**

**20mm**

**MFE-821-30**


**30mm**

**MFE-821-40**

**40mm**

**MFE-821-60**

**60mm**




**Cottle Speculum**



- 85mm length

**Code**      **Blunt Prongs**

<b>MFE-822-S</b>	<b>19cm/7"</b>	<b>Straight</b>
<b>MFE-822-R</b>	<b>19cm/7"</b>	<b>Cranked Right</b>
<b>MFE-822-L</b>	<b>19cm/7"</b>	<b>Cranked Left</b>



**Kilner Alae Retractor**

- 2 prong blunt ball tips
- Knife guard

**Code**      **Length**

<b>MFE-823-R</b>	<b>125mm x</b>	<b>Cranked Right</b>
<b>MFE-823-L</b>	<b>85mm</b>	<b>Cranked Left</b>



**Joseph Nasal Saw**

- 19.5cm/7½" length
- Bayonet shaped

**Code**      **Length**

<b>MFE-824-03</b>	<b>3mm</b>
<b>MFE-824-04</b>	<b>4mm</b>
<b>MFE-824-05</b>	<b>5mm</b>



**Ballenger Swivel Knife**

- Double Action

**Code**      **Length**

<b>MFE-825-19</b>	<b>19cm/7½"</b>
-------------------	-----------------




**Jansen Middleton Septum Forceps**



- 23cm/9" length

Code	Blade
MFE-826-U	Universal
MFE-826-L	Left
MFE-826-R	Right



**Walsham Septum Straightening Forceps**

Code	Length
MFE-827-23	23cm/9"



**Asch Septum Straightening Forceps**



- Head size 22mm Ø
- Weight 165g/6oz

Code	Length
MFE-828-19	19cm/7½"



**Toffee Mallet**

- Head size 20mm Ø
- Weight 145g/5oz
- Tufnol handle

Code	Length
MFE-829-17	17cm/6¾"



**Ultra Fine Mallet**

- 18cm/7" length

Code	Width
MFE-830-04	4mm
MFE-830-05	5mm
MFE-830-06	6mm



**Read Nasal Gouge**

- Round cutting edge 6mm width
- 16cm/6¼" length
- Bayonet shaped
- Notched cutting edge 6mm width

Code	Length
MFE-831-R	Round 6mm
MFE-831-N	Notched 6mm



**Killian Septum Gouge**

- 18cm/7¼" length
- Graduated

Code	Width
MFE-832-04	4mm
MFE-832-06	6mm
MFE-832-08	8mm



**Obwegeser Nasal Septum Osteotome**



- 11mm wide
- Bevelled at one corner

Code	Length
MFE-834-18	18cm/7"



**Matthew Osteotome**

- 3mm wide
- Aluminium handle

Code	Length
MFE-835-18	18.5cm/7¼"



**Gubisch Osteotome**

- 18cm/7" Length

Code	Size
MFE-836-08	8mm
MFE-836-10	10mm
MFE-836-12	12mm
MFE-836-14	14mm
MFE-836-16	16mm



**Rubin Osteotome**

- 14cm/5½" length

Code	Size
MFE-837-02	2mm
MFE-837-05	5mm
MFE-837-07	7mm
MFE-837-09	9mm
MFE-837-11	11mm
MFE-837-13	13mm



**French (Gillies) Osteotome**

- 16cm/6¼" length

Code	Size
MFE-838-08	8mm
MFE-838-10	10mm
MFE-838-12	12mm
MFE-838-14	14mm
MFE-838-16	16mm



**Cinelli Double Guarded Osteotome**



- 18cm/7" length

Code	Size
MFE-839-02	2mm
MFE-839-04	4mm
MFE-839-05	5mm
MFE-839-06	6mm
MFE-839-08	8mm
MFE-839-10	10mm
MFE-839-12	12mm



**Jayes Osteotome**

- 18.5cm/7½" length

Code	Size
MFE-840-02	2mm
MFE-840-04	4mm
MFE-840-05	5mm
MFE-840-06	6mm
MFE-840-08	8mm
MFE-840-10	10mm
MFE-840-12	12mm



**Cottle Chisel**

- 15cm/6" length

Code	Size
MFE-841-02	2mm
MFE-841-04	4mm
MFE-841-05	5mm
MFE-841-06	6mm
MFE-841-08	8mm
MFE-841-10	10mm
MFE-841-12	12mm



**McIndoe Chisel**

- 15cm/6" length

Code	Size
MFE-842-02	2mm
MFE-842-04	4mm
MFE-842-05	5mm
MFE-842-06	6mm
MFE-842-08	8mm
MFE-842-10	10mm
MFE-842-12	12mm



- Set of 5

Code	Sizes
MFE-845-02	2mm to 10mm

**Sheehan Osteotomes**



- 14.5cm/5<sup>3</sup>/<sub>4</sub>" length
- 5mm wide

**Code**

<b>MFE-846</b>	Without guard
<b>MFE-846-L</b>	Guard on left
<b>MFE-846-R</b>	Guard on right



**Braithwaite Chisel**

- 16.5cm/6<sup>1</sup>/<sub>2</sub>" length
- Guard on Left

**Code**

**Size**

<b>MFE-847-04</b>	<b>4mm</b>
<b>MFE-847-05</b>	<b>5mm</b>
<b>MFE-847-06</b>	<b>6mm</b>
<b>MFE-847-07</b>	<b>7mm</b>



**Foman Chisel with Guard**

- 18cm/7" length

**Code**

<b>MFE-848-S</b>	Straight
<b>MFE-848-L</b>	Curved left
<b>MFE-848-R</b>	Curved right

**Silver Chisel with Guard (Frenchay)**

- 20cm/8" length

**Code**

**Size**


<b>MFE-849-S</b>	Straight
<b>MFE-849-L</b>	Curved left
<b>MFE-849-R</b>	Curved right



**Nievert Chisel with Guard**



- 21cm/8¼” length
- Tungsten Carbide surface




**Code**

<b>MFE-850-C</b>	<b>Course</b>
<b>MFE-850-M</b>	<b>Medium</b>
<b>MFE-850-F</b>	<b>Fine</b>

**Fomon Double Ended Nasal Rasp**

- 18cm/7” length
- Tungsten Carbide
- surface Backward cutting



**Code**

**Size**

<b>MFE-851-C</b>	<b>Course</b>
<b>MFE-851-M</b>	<b>Medium</b>
<b>MFE-851-F</b>	<b>Fine</b>

**Foman Chisel with Guard**



**Code**

**Length**

**MFE-852-18**

**18.5cm/7½"**



**Maltz Nasal Rasp**

- Tungsten Carbide surface

**Code**

**Length**

**MFE-853-18**

**18.5cm/7½"**



**Maltz Nasal Rasp**

- Diamond surface

**Code**

**Tip**

**Length**

**MFE-854-C**

**Curved**

**17.5cm/6¾"**

**MFE-854-S**

**Straight**



**Peet Nasal Rasp**

- Diamond surface

**Code**

**Length**

**MFE-855-22**

**22cm/8¾"**



**Aufrecht (Peet) Nasal Rasp**



- Tungsten Carbide surface

Code	Length
<b>MFE-856-21</b>	<b>21cm/8<sup>1</sup>/<sub>4</sub>"</b>



**Glabella Nasal Rasp**

- Diamond surface


Code	Length
<b>MFE-857-21</b>	<b>21cm/8<sup>1</sup>/<sub>4</sub>"</b>



**Fomon Double Ended Nasal Rasp**

- Tungsten Carbide surface

Code	Length
<b>MFE-858-18</b>	<b>18cm/7"</b>



**Miller Nasal Rasp**



For endoscopic forehead face lift surgery with ergonomic handles.

- Straight 10mm wide
- For initial dissection through forehead incisions

Code	Length
<b>MFE-859-22</b>	22.5cm/8 <sup>3</sup> / <sub>4</sub> "



**Frontotemporal Dissector**

- Curved 7mm wide
- For initial dissection

Code	Length
<b>MFE-860-24</b>	24cm/9 <sup>1</sup> / <sub>2</sub> "



**Frontotemporal Dissector**

- Curved 9mm wide
- For a posterior scalp dissection

Code	Length
<b>MFE-861-25</b>	25cm/9 <sup>3</sup> / <sub>4</sub> "



**Periosteal Dissector**

- S-Shaped 7mm wide
- Used in the dissection of the naso-glabbellar area of the forehead

Code	Length
<b>MFE-862-26</b>	26cm/10 <sup>1</sup> / <sub>2</sub> "



**Frontoglabbellar Dissector**



For endoscopic forehead face lift surgery with ergonomic handles.

- Straight 10mm wide
- Separates the temporal fascia



Code	Length
<b>MFE-863-23</b>	<b>23.5cm/9<sup>1</sup>/<sub>4</sub>"</b>

**Temporal Line Dissector**

- Straight 5mm wide
- For use in the lower eyelid approach



Code	Length
<b>MFE-864-22</b>	<b>22.5cm/8<sup>3</sup>/<sub>4</sub>"</b>

**Zygomatic Arch Elevator**

- Curved 12mm wide



Code	Length
<b>MFE-865-25</b>	<b>25.5cm/10"</b>

**Scalp Elevator**

- Half curved with angled tip 12mm wide



Code	Length
<b>MFE-866-24</b>	<b>24cm/9<sup>1</sup>/<sub>2</sub>"</b>

**Scalp Elevator**



For endoscopic forehead face lift surgery with ergonomic handles.

- Half curved with angled tip  
12mm wide



Code	Length
MFE-867-24	24cm/9½”

**Scalp Elevator**

- Curved 6mm wide
- For dissection of the deep temporal fascia



Code	Length
MFE-868-22	22.5cm/8¾”

**Flap Dissector**

- S-Shaped 8mm wide
- For dissection along the orbital rim



Code	Length
MFE-869-24	24cm/9½”

**Arcus Marginalis Dissector**

- Half curved 7mm wide



Code	Length
MFE-870-24	24cm/9½”

**Endoforehead Elevator**



For endoscopic forehead face lift surgery with ergonomic handles.

- Curved tip 7mm wide

Code	Length
<b>MFE-871-23</b>	<b>23.5cm/9¼"</b>



**Transoral Dissector**

- Half curved

Code	Length
<b>MFE-872-24</b>	<b>24cm/9½"</b>



**Orbital Rim Dissector**

Code	Length
<b>MFE-873-24</b>	<b>24cm/9½"</b>



**Nerve Hook - Curved Right**

Code	Length
<b>MFE-874-24</b>	<b>24cm/9½"</b>



**Nerve Hook - Curved Left**



For endoscopic forehead face lift surgery with ergonomic handles.

**Code**

**Length**

**MFE-875-17**

**17cm/9<sup>3</sup>/<sub>4</sub>"**



**Nerve Retractor - Curved Left**

**Code**

**Length**

**MFE-876-17**

**17cm/9<sup>3</sup>/<sub>4</sub>"**



**Nerve Retractor - Curved Right**

- Slightly curved Round edge tip
- 8mm wide

**Code**

**Length**

**MFE-877-26**

**26cm/10<sup>1</sup>/<sub>4</sub>"**



**Suction Elevator**

- Half curved Round edge tip
- 8mm wide

**Code**

**Length**

**MFE-878-26**

**26cm/10<sup>1</sup>/<sub>4</sub>"**



**Suction Elevator**



For endoscopic forehead face lift surgery with ergonomic handles.

- Slightly curved

Code	Length
MFE-879-18	18cm/7"



**Mini Forehead Elevator**

- Half curved

Code	Length
MFE-880-18	18cm/7"



**Mini Forehead Elevator**

- Fully Curved

Code	Length
MFE-881-18	18cm/7"



**Mini Forehead Dissector**



- 18cm/7" length
- Straight

Code      Width

**MFE-882-04**      4mm

**MFE-882-06**      6mm

**MFE-882-08**      8mm

- 17.5cm/6¾" length

- Curved

Code      Width

**MFE-883-04**      4mm

**MFE-883-06**      6mm

**MFE-883-08**      8mm



**Epkers Osteotome**

- 16cm/6¼" length

- Straight

Code      Width

**MFE-884-02**      2.5mm

**MFE-884-03**      3.5mm

**MFE-884-05**      5mm

**MFE-884-07**      7mm

**MFE-884-10**      10mm

**MFE-884-15**      15mm

**MFE-884-18**      18mm



**Tessier Multi-Purpose Osteotome**

- 16cm/6¼" length

- Short blade with rounded edges

- Straight

Code      Width

**MFE-885-10**      10mm

**MFE-885-15**      15mm

- Curved

Code      Width

**MFE-886-10**      10mm

**MFE-886-15**      15mm



**Tessier Bone Graft Osteotome**



- 23.5cm/9" length
- Curved

Code	Size
MFE-887-08	8mm
MFE-887-10	10mm
MFE-887-11	11mm



**Obwegeser Pterygoid Osteotome**

- 22cm/8<sup>3</sup>/<sub>4</sub>" length
- Thin blade

Code	Size
MFE-888-08	8mm
MFE-888-12	12mm
MFE-888-16	16mm



**Obwegeser Splitting Osteotome**

- 22cm/8<sup>3</sup>/<sub>4</sub>" length
- Wedge profile

Code	Size
MFE-889-08	8mm
MFE-889-12	12mm
MFE-889-16	16mm



**Obwegeser-Ramus Wedge Osteotome**

- 22cm/8<sup>3</sup>/<sub>4</sub>" length
- Thicker tapered blade

Code	Size
MFE-890-07	7mm
MFE-890-10	10mm



**Ramus Split Osteotome Straight**



**Disimpaction Forceps**

Nasal blade

Oral blade


Oral blade

Nasal blade



- Tessier-Rowe Left
- For Children

Code	Size
MFE-891-23	23cm/9"




- Tessier-Rowe Right
- For Children

Code	Size
MFE-893-23	23cm/9"




- Tessier-Rowe Left
- For Adult

Code	Size
MFE-892-23	23cm/9"



- Tessier-Rowe Right
- For Adult

Code	Size
MFE-894-23	23cm/9"



Nasal blade

Oral blade

**Tessier-Rowe Left Side**


Oral blade

Nasal blade



These instruments are used for maxillary advancement in Le Fort I, II, or III procedures. A Large curvature of the blade prevents damage of the teeth. The oral blades are asymmetrical to adapt to the curvature of the palate which is often high and concave.

**Tessier-Rowe Right Side**






**Disimpaction Forceps**

Nasal blade

Oral blade

Oral blade


Nasal blade



- Tessier-Rowe Left
- For Children

Code      Size

MFE-895-23    23cm/9"



- Tessier-Rowe Right
- For Children

Code      Size


MFE-897-23    23cm/9"



- Tessier-Rowe Left
- For Adult

Code      Size

MFE-896-23    23cm/9"



- Tessier-Rowe Right
- For Adult

Code      Size

MFE-898-23    23cm/9"

Nasal blade

Oral blade

Tessier-Rowe Left Side

Oral blade

Nasal blade

40mm

These instruments are used for maxillary advancement in Le Fort I, II, or III prodecures. A Large curvature of the blade prevents damage of the teeth. The oral blades are asymmetrical to adapt to the curvature of the palate which is often high and concave.

Tessier-Rowe Right Side


40mm



## Smith Ramus Separator

- 12mm width

Code	Size
MFE-899-23	23cm/9"



The superior ramus separator in conjunction with the sagittal split separator allows controlled separation of the superior ramus osteotomy cut. It's unique blade action exerts force over a broad surface area for maximum control. This instrument is also useful in "down facturing" the maxilla by placement at the zygomatic buttress osteotomy cut.




## Obwegeser Fragment Adapting Forceps Left Side

**Code**

**Size**

**MFE-900-19**

**19.5cm/7<sup>3</sup>/<sub>4</sub>"**



## Obwegeser Fragment Adapting Forceps Right Side

**Code**

**Size**

**MFE-901-19**


**19.5cm/7<sup>3</sup>/<sub>4</sub>"**





- Tips 0.2mm Ø

Code	Size
<b>MFE-902-11</b>	11cm/4 <sup>1</sup> / <sub>4</sub> "



**Vessel Dilator Angled**

- Tips 0.3mm Ø



Code	Size
<b>MFE-903-11</b>	11cm/4 <sup>1</sup> / <sub>4</sub> "



**Vessel Dilator Angled**

- Tips 0.2mm Ø


Code	Size
<b>MFE-904-13</b>	13.5cm/5 <sup>3</sup> / <sub>8</sub> "



**Vessel Dilator**

- Curved
- Tips 0.2mm Ø

Code	Size
<b>MFE-905-13</b>	13.5cm/5 <sup>3</sup> / <sub>8</sub> "




**Vessel Dilator**



- Tips 0.3mm Ø

Code	Size
<b>MFE-906-18</b>	<b>18cm/7"</b>




**Vessel Dilator Long**

- Curved
- Tips 0.2mm Ø

Code	Size
<b>MFE-907-18</b>	<b>18cm/7"</b>



**Vessel Dilator Long**




Mast Pak Surgical Corp. Vessel Dilators have Mirror polished parallel tips for interluminal vessel dilation. The entire area is opened to the same diameter and the pressure is distributed evenly. This ensures that unnecessary damage to the vessel is avoided.



- Straight
- Tips 0.3mm Ø Balanced forceps


Code	Size
MFE-908-15	15cm/6"



**Round Bodied Vessel Dilator**

- Curved
- Tips 0.2mm Ø Balanced forceps


Code	Size
MFE-909-15	15cm/6"



**Round Bodied Vessel Dilator**

- Straight
- Tips 0.3mm Ø Balanced forceps

Code	Size
MFE-910-18	18cm/7"



**Round Bodied Vessel Dilator**

- Curved
- Tips 0.3mm Ø Balanced forceps

Code	Size
MFE-911-18	18cm/7"



**Round Bodied Vessel Dilator**



- Tips 0.3mm

Code	Size
MFE-912-11	11cm/4 <sup>1</sup> / <sub>4</sub> "

**Jewellers Forceps**



- Tips 0.3mm

Code	Size
MFE-913-15	15cm/6"

**Jewellers Forceps Long**



- Tips 0.3mm

Code	Size
MFE-914-12	12cm/4 <sup>3</sup> / <sub>4</sub> "

**Jewellers Forceps**



- Tips 0.3mm

Code	Size
MFE-915-13	13.5cm/5 <sup>3</sup> / <sub>8</sub> "

**Jewellers Forceps**



- Tips 0.3mm

Code	Size
MFE-916-12	12cm/4 <sup>3</sup> / <sub>4</sub> "

**Pierse Forceps**

Pierse Forceps have a micro hole in the tips which securely grip tissue with minimal trauma



- 1 x 2 micro teeth

- Tips 0.3mm

Code	Size
MFE-917-12	12cm/4 <sup>3</sup> / <sub>4</sub> "

**Micro Adson Fixation Forceps**





- Straight, with platform
- Tips 0.3mm Ø Balanced forceps

**Code**                      **Size**

**MFE-918-15**    **15cm/6"**



**Round Bodied Forceps**

- Curved, with platform
- Tips 0.3mm Ø Balanced forceps

**Code**                      **Size**

**MFE-919-15**    **15cm/6"**



**Round Bodied Forceps**

- Angled
- Tips 0.2mm Ø Balanced forceps

**Code**                      **Size**

**MFE-920-15**    **15cm/6"**



**Round Bodied Angled Forceps**

- Straight, with platform
- Tips 0.5mm Ø Balanced forceps

**Code**                      **Size**

**MFE-921-18**    **18cm/7"**



**Round Bodied Forceps**

- Curved, with platform
- Tips 0.5mm Ø Balanced forceps

**Code**                      **Size**

**MFE-922-18**    **18cm/7"**



**Round Bodied Forceps**

- Angled
- Tips 0.2mm Ø Balanced forceps

**Code**                      **Size**

**MFE-923-18**    **18cm/7"**




**Round Bodied Angled Forceps**



- Blades 12mm long
- Blades 14mm long
- Adventitia scissors with pointed tips


Code	Size
MFE-924-15	15cm/6"
MFE-924-18	18cm/7"



**Straight Scissors Round Handle**

- Blades 12mm long
- Blades 14mm long
- Dissecting scissors with rounded tips

Code	Size
MFE-925-15	15cm/6"
MFE-925-18	18cm/7"



**Curved Scissors Round Handle**

- Blades 15mm long Curved blades

Code	Size
MFE-926-15	15cm/6"
MFE-926-18	18cm/7"



**Blondeel Scissors Round Handle**

These Scissors were designed by professor Phillip Blondeel for use in breast reconstruction utilising a perforator ap. A combination of ring handle design for preparation of ne structures, the round handle micro instrument handle allows the user to perform ne micro surgery.



- Round bodied 0.4mm jaws

**Without Lock**

**Curved**

**Straight**

13cm/5"

MFE-927-13

MFE-928-13

15cm/6"

MFE-927-15

MFE-928-15

18cm/7"

MFE-927-18

MFE-928-18



**Barraquer Micro Needle Holder**

- Flat handled 0.4mm jaws

**Without Lock**

**Curved**

14cm/5½"

MFE-929-14

18cm/7"

MFE-929-18



**Castroviejo Micro Needle Holder**

- Round bodied 0.4mm jaws

**Without Lock**

**Curved**

15cm/6"

MFE-930-15

18cm/7"

MFE-930-18




**Double Action Micro Needle Holder**



- Tips 0.1mm Ø


Code	Size
MFE-931-11	11cm/4 <sup>1</sup> / <sub>4</sub> "



**Ultra Fine Jewellers Forceps**

- Tips 0.1mm Ø


Code	Size
MFE-932-13	13.5cm/5 <sup>3</sup> / <sub>8</sub> "



**Ultra Fine Jewellers Forceps**

- Straight
- Tips 0.1mm Ø Balanced forceps


Code	Size
MFE-933-15	15cm/6"



**Ultra Fine Round Bodied Forceps**

- Straight
- Tips 0.1mm Ø Balanced forceps

Code	Size
MFE-934-18	18cm/7"



**Ultra Fine Round Bodied Forceps**



- Tips 0.1mm Ø

Code	Size
MFE-935-11	11cm/4 <sup>1</sup> / <sub>4</sub> "



**Ultra Fine Vessel Dilator Angled**

- Curved tips 0.2mm Ø


Code	Size
MFE-936-13	13.5cm/5 <sup>3</sup> / <sub>8</sub> "



**Ultra Fine Vessel Dilator Angled**

- Straight tips 0.1mm Ø Balanced forceps


Code	Size
MFE-937-15	15cm/6"



**Ultra Fine Round Bodied Vessel Dilator**

- Curved tips 0.2mm Ø Balanced forceps

Code	Size
MFE-938-15	15cm/6"



**Ultra Fine Round Bodied Vessel Dilator**



- Blades 14mm long
- Adventitia scissors with pointed tips

**Code**

**Size**

**MFE-939-15**

**15cm/6"**



**Ultra Fine Straight Scissors Flat Handle**



- Blades 12mm long
- Dissecting scissors with rounded tips

**Code**

**Size**

**MFE-940-15**

**15cm/6"**



**Ultra Fine Curved Scissors Flat Handle**



- Blades 14mm long
- Dissecting scissors with pointed tips

**Code**

**Size**

**MFE-941-15**

**15cm/6"**

**MFE-941-18**

**18cm/7"**



**Ultra Fine Straight Scissors Round Handle**



- Blades 12mm long
- Dissecting scissors with rounded tips

**Code**

**Size**

**MFE-942-15**

**15cm/6"**

**MFE-942-18**

**18cm/7"**



**Ultra Fine Curved Scissors Round Handle**





• Flat handled 0.2mm Jaws

Code	Size
MFE-943-14	14cm/5½"



**Ultra Fine Curved Micro Needle Holder**

• Round bodied 0.2mm Jaws

Code	Size
MFE-944-15	15cm/6"
MFE-944-18	18cm/7"



**Ultra Fine Curved Micro Needle Holder**



**INSTRUMENTS CONTAINERS**



**SINGLE TRAY CONTAINERS**

Internal Dimensions			Code
A	B	C	
200 x 140 x 50mm			<b>MFE-945-200</b>
280 x 195 x 50mm			<b>MFE-945-280</b>
265 x 255 x 50mm			<b>MFE-945-265</b>
380 x 260 x 50mm			<b>MFE-945-380</b>

**Standard Container illustrated**  
One of our popular cases for  
Micro Instruments



**DOUBLE TRAY CONTAINERS**

Internal Dimensions			Code
A	B	C	
280 x 195 x 65mm			<b>MFE-946-280</b>
265 x 255 x 65mm			<b>MFE-946-265</b>
380 x 260 x 65mm			<b>MFE-946-380</b>

**Double Tray Container illustrated**  
Two removable perforated  
trays



**Acland Clamp Container**

With individual compartments for  
various clamp sizes

**Code**  
**MFE-947**

## More Information



Visit Our Website  
[www.rmfelahi.com](http://www.rmfelahi.com)



Our Location  
**NEAR ROUND ABOUT  
S.I.E, SIALKOT-  
PAKISTAN**



[export@mfelahi.com](mailto:export@mfelahi.com)  
[zakriamfe@gmail.com](mailto:zakriamfe@gmail.com)



+92-333-8664664  
+92-336-7322733  
+92-333-8384444

